


Product Datasheet

Rat Endostatin ELISA Kit (orb381124)

Description	Rat Endostatin ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: H1591-K1774
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:5, concentration in serum and plasma is around 10 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	18933 MW
Uniprot ID	Q9WUW5
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Rat
Endostatin
PicoKine
ELISA Kit
standa...