


Product Datasheet

Mouse TPOR/MPL/Thrombopoietin R ELISA Kit (orb381129)

Description	Mouse TPOR/MPL/Thrombopoietin R ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q26-W482
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	69788 MW
Uniprot ID	Q08351
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Mouse TPOR/MPL PicoKine ELISA Kit standa...