


Product Datasheet

Mouse FCRN/FCGRT ELISA Kit (orb381131)

Description	Mouse FCRN/FCGRT ELISA Kit
Reactivity	Mouse
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S22-S297
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is 100-1000p/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
MW	40093 MW
Uniprot ID	Q61559
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
FCRN/FCGRT
PicoKine
ELISA Kit
stan...