


Product Datasheet

Human EDIL3 ELISA Kit (orb381149)

Description	Human EDIL3 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D24-E480
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 150 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	53765 MW
Uniprot ID	043854
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
EDIL3
PicoKine
ELISA Kit
standard ...