


## Product Datasheet

### Rat GDF5/Bmp 14 ELISA Kit (orb381153)

<b>Description</b>	Rat GDF5/Bmp 14 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	31.2 pg/ml - 2,000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: A376-R495
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml
<b>Sample Types</b>	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
<b>MW</b>	7491 MW
<b>Uniprot ID</b>	<b>Q8K4X5</b>
<b>Sensitivity</b>	10 pg/ml
<b>Expiration Date</b>	Please enquire.



Rat GDF5  
PicoKine  
ELISA Kit  
standard  
cur...