


Product Datasheet

Human Siglec-7/CD328 ELISA Kit (orb381155)

Description	Human Siglec-7/CD328 ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q19-L353
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	51143 MW
Uniprot ID	Q9Y286
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
Siglec-
7/CD328
PicoKine
ELISA Kit ...