


Product Datasheet

Mouse ALK-1/ACVRL1 ELISA Kit (orb381161)

Description	Mouse ALK-1/ACVRL1 ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D23-P119
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	56519 MW
Uniprot ID	Q61288
Sensitivity	15 pg/ml
Expiration Date	Please enquire.



Mouse ALK-1/ACVRL1 PicoKine ELISA Kit st...