


Product Datasheet

Human TREML1/TIt 1 ELISA Kit (orb381163)

Description	Human TREML1/TIt 1 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q16-P162
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is around 30 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	32679 MW
Uniprot ID	Q86YW5
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
TREML1
PicoKine
ELISA Kit
standard...