



Product Datasheet

USO1 antibody (orb382113)


Description	Rabbit polyclonal antibody to USO1.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IH, IP, WB
Immunogen	KLH-conjugated synthetic peptide of human p115
Target	USO1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O60763
Entrez	8615
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IP: 1:20-50
Expiration Date	12 months from date of receipt.



Western blot analysis of p115 expression...



Immunohistochemical analysis of p115 sta...



Immunofluorescent analysis of p115 stain...