
Product Datasheet

B7-H6 antibody (orb382134)

Description

B7-H6 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human B7-H6. Please contact us for the exact immunogen sequence. The peptide is available as orb13013.

Target

B7-H6

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW


51 kDa

Uniprot ID
[Q68D85](#)
NCBI
[001189368](#), [001202439](#)
Entrez
[374383](#)
Dilution Range


IHC-P: 1:50-300, WB: 1:200-1000

Expiration Date


12 months from date of receipt.



IHC-P image of guinea pig brain tissue u...



Immunohistochemical staining of paraffin...



Immunohistochemical staining of human br...