



Product Datasheet

RANTES antibody (orb382138)


Description	RANTES antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein of human RANTES
Target	RANTES
Form/Appearance	PBS with 0.05% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified by protein G.
MW	30 kDa
Uniprot ID	P30882, P13501, P50231
Entrez	51621
Dilution Range	WB: 1:100-2000, IHC-P: 1:100-300
Expiration Date	12 months from date of receipt.



IHC-P staining of using RANTES antibody ...



Immunohistochemical staining of using RA...



Immunohistochemical staining of paraffin...