

Product Datasheet

CREB1 peptide (orb382236)

Description	The peptide is used to block anti-CREB1 antibody reactivity
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	CREB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382519 (CREB1 antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q01147 , P15337 , P16220
Entrez	1385 , 12912 , 81646
Expiration Date	12 months from date of receipt.