

Product Datasheet

MRPS18C peptide (orb382264)

Description	The peptide is used to block anti-MRPS18C antibody reactivity
Reactivity	Human, Primate
Tested Applications	Blocking
Target	MRPS18C
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382547 (MRPS18C antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Y3D5
Entrez	51023
Expiration Date	12 months from date of receipt.