

Product Datasheet

KRT2 peptide (orb382315)

Catalog Number	orb382315
Category	Proteins
Description	The peptide is used to block Anti-Cytokeratin-pan Antibody (orb382594) reactivity.
Reactivity	Bovine, Canine, Human, Mouse, Rat
Form/Appearance	Lyophilized powder
Purity	>85%
UniProt ID	Q86Y46, P04259, P48668, P12035, P35908, O95678, P08729, Q01546, P13647, Q14CN4, Q9NSB2, Q922U2, Q6IG02, Q5XKE5, Q3SY84, P02538, Q3UV17, P05787, Q6P6Q2, Q3TTY5, Q7RTS7
Tested applications	Blocking
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Source	Synthetic
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Entrez	77055, 3852, 319101, 369017, 3854, 338785, 121391, 16681, 3850, 51350, 140807, 406228, 3890, 286887, 112802, 3853, 110308, 3849, 3856, 9119

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)