

Product Datasheet

TUBA1A peptide (orb382316)

Description	The peptide is used to block anti-TUBA1A antibody reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TUBA1A; TUBA1B; TUBA1C; TUBA3C; TUBA4A; TUBA8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb382595 (TUBA1A antibody)Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9BQE3 , Q6AYZ1 , P68370 , Q9NY65 , P05213 , Q6P9V9 , Q71U36 , P68363 , Q13748 , P68373 , P68366 , P68369
Entrez	7846 , 7277 , 10376 , 22146 , 64158 , 22142 , 84790 , 7278 , 113457 , 100909441 , 51807 , 1962 , 500929 , 22143
Expiration Date	12 months from date of receipt.