


Product Datasheet

CD44 protein (orb382874)

Description	CD44 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	CD44
Form/Appearance	50mM Tris,50mM NaCl,2M Urea,1mM DDT,pH8.0
Concentration	0.5 mg/ml
Storage	Store at -20°C for at least 1 year. After reconstitution, protein solution is stable at -20°C for 6 months to 12 months. Keep at -80°C for long term storage. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	N-Trx-His6
Purity	> 85%
MW	33.7 kDa
NCBI	004372
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of CD44 protein