




Product Datasheet

FGFR4 antibody (orb38290)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	FGFR4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87954
Uniprot ID	P22455
NCBI	NP_998812.1 , NP_002002.3 , NP_075252.2
Dilution Range	WB: 1:1000, FACS: 1:10-50, IF/ICC: 1:200, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of FGFR4 antibody ...



Flow cytometric analysis of FGFR4 antibo...