

Product Datasheet

Human FCER1G protein (orb382950)

Description	Recombinant human FCER1G protein.
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	E.coli and Yeast N-terminal GST-tagged Partial
Protein Sequence	LGEPQLCYILDAILFLYGIVLTLTYCRLKIQRKAAITSYEKSDGVYTGLSTRNQETY ETLKHEKPPQ
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	34.8 kDa
Uniprot ID	P30273
Expiration Date	6 months from date of receipt.