

Product Datasheet


Cynomolgus M-CSF R / CSF1R / CD115 Protein (orb383516)

Description


Cynomolgus M-CSF R, His Tag is expressed from



Reactivity	Pro ...
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	M-CSF R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.1 kDa. The protein migrates as 70-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Cynomolgus M-CSF R, His Tag (orb383516) is expressed from human 293 cells (HEK293). It contains AA Ile 99 - Pro 596 (Accession # G7P8P1).
MW	57.1 kDa
Uniprot ID	G7P8P1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Monkey M-CSF protei...



Immobilized Cynomolgus M-CSF R, His Tag ...