


Product Datasheet

AMBRA1 Antibody (orb384088)

Description	Rabbit Anti-Human AMBRA1 Polyclonal
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide of Human AMBRA1 (aa. 200-300)
Target	AMBRA1
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-644 was sufficient for detection of AMBRA1 in 15 µg of Rat Brain Lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	142 kDa
Uniprot ID	Q9C0C7
NCBI	NP_001254711.1
Entrez	55626
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.



Western Blot analysis of rat brain cell

...