

Product Datasheet

ALG12 antibody (HRP) (orb38460)

Description	ALG12 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Dol-P-Man:Man (7)GlcNAc (2)-PP-Dol alpha-1,6-mannosyltransferase protein (367-422AA)
Target	ALG12
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BV10
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.