


Product Datasheet

Human IGFBP4 ELISA Kit (orb385306)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------|
| Description | Sandwich ELISA for quantitative detection of human IGFBP4 concentrations. |
| Reactivity | Human |
| Range | 500 pg/ml - 32000 pg/ml |
| Tested Applications | ELISA |
| Target | IGFBP4 |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | P22692 |
| Sensitivity | 20 pg/ml |
| Entrez | 3487 |
| Expiration Date | Please enquire. |



Representative Standard Curve