




Product Datasheet

CCK4 antibody (orb38531)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PTK7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	118392
Uniprot ID	Q13308
NCBI	NP_690619.1 , NP_002812.2 , NP_690620.1 , NP_690621.1 , NP_001257327.1
Dilution Range	FACS: 1:10-50, IHC-P: 1:50-100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of CCK4 antibody...



Immunohistochemical analysis of formalin...