


Product Datasheet

Human NECTIN4 ELISA Kit (orb385321)

Description	Sandwich ELISA for quantitative detection of human Nectin 4 concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	NECTIN4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates, Urine
Uniprot ID	Q96NY8
Sensitivity	10 pg/ml
Entrez	81607
Expiration Date	Please enquire.



Representative
Standard
Curve