

Product Datasheet

Human PEDF ELISA Kit (orb385340)

Description Sandwich ELISA for quantitative detection of human PEDF concentrations.

Reactivity Human

Range 780 pg/ml - 50000 pg/ml

Tested Applications ELISA

Target SERPINF1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P36955**

Sensitivity 20 pg/ml

Entrez **5176**

Expiration Date Please enquire.



Representative
Standard
Curve