


Product Datasheet

Human SHBG ELISA Kit (orb385342)

Description	Sandwich ELISA for quantitative detection of human SHBG concentrations.
Reactivity	Human
Range	0.031 nmol/L - 2 nmol/L
Tested Applications	ELISA
Target	SHBG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P04278
Sensitivity	0.001 nmol/L
Entrez	6462
Expiration Date	Please enquire.



Representative Standard Curve