


Product Datasheet

Mouse ACVRL1 ELISA Kit (orb385358)

Description	Sandwich ELISA for quantitative detection of mouse ALK1 concentrations.
Reactivity	Mouse
Range	46.9 pg/ml - 3000 pg/ml
Tested Applications	ELISA
Target	ACVRL1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q61288
Sensitivity	15 pg/ml
Entrez	11482
Expiration Date	Please enquire.



Representative
Standard
Curve