


Product Datasheet

Mouse CXCL9 ELISA Kit (orb385371)

Description	Sandwich ELISA for quantitative detection of mouse CXCL9 concentrations.
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	CXCL9
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P18340
Sensitivity	10 pg/ml
Entrez	17329
Expiration Date	Please enquire.



Representative
Standard
Curve