


Product Datasheet

Mouse LEP ELISA Kit (orb385387)

Description	Sandwich ELISA for quantitative detection of mouse Leptin concentrations.
Reactivity	Mouse
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	LEP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P41160
Sensitivity	60 pg/ml
Entrez	16846
Expiration Date	Please enquire.



Representative
Standard
Curve