


Product Datasheet

Rat CXCL5 ELISA Kit (orb385416)

Description	Sandwich ELISA for quantitative detection of rat CXCL5 concentrations.
Reactivity	Rat
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	CXCL5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P97885
Sensitivity	1 pg/ml
Entrez	60665
Expiration Date	Please enquire.



Representative Standard Curve