



## Product Datasheet

### SLC22A8 antibody (orb385519)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SLC22A8                                                                                                      |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | WB                                                                                                                                         |
| <b>Immunogen</b>           | A synthetic peptide of human SLC22A8                                                                                                       |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.                                                                                 |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks.<br>For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                 |
| <b>Purity</b>              | Affinity purification                                                                                                                      |
| <b>MW</b>                  | 60kDa                                                                                                                                      |
| <b>Uniprot ID</b>          | <a href="#">Q8TCC7</a>                                                                                                                     |
| <b>Entrez</b>              | <a href="#">9376</a>                                                                                                                       |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:25-100                                                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |



Western blot analysis of extracts of var...



Immunohistochemical staining of human br...