



Product Datasheet

TNF alpha Antibody (orb385905)


Description	Tumor Necrosis Factor Alpha (TNF alpha) is a cytokine that binds to TNFRSF1A/TNFR1 and TNFRSF1B/TNFB...
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FACS, IF, IHC-P
Immunogen	Full length human recombinant protein was used as the immunogen for the TNF alpha antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the TNF alpha antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Titering of the TNF alpha antibody may be required for optimal performance.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgM, kappa
Clonality	Monoclonal
Purity	PEG precipitation
Uniprot ID	P01375
Hazard Information	This TNF alpha antibody is available for research use only.
Dilution Range	Immunofluorescence: 1-2ug/ml, Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 2-4ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Immunofluorescent staining of PFA-fixed ...



IHC testing of FFPE human histiocytoma a...



IHC testing of FFPE rat stomach and TNF ...