



Product Datasheet

Aurora B Antibody / AURKB (orb385982)

Description	Recognizes a protein of 39-45kDa, which is identified as Aurora B. The serine/threonine protein kina...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Amino acids 89-251 were used as the immunogen for this Aurora B antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide
Storage	Store the Aurora B antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Aurora B antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity
Uniprot ID	Q96GD4
Hazard Information	This Aurora B antibody is available for research use only.
Entrez	9212
Dilution Range	Western blot: 0.5-1ug/ml
Expiration Date	12 months from date of receipt.



Western blot testing of 1) a partial rec...



SDS-PAGE Analysis of Purified, BSA-Free ...