

---

## Product Datasheet

**TNFSF15 Antibody / VEGI / TL1A (orb386086)**

**Description**

TNF ligand-related molecule 1 is an anti-angiogenic

**Species/Host**

Mouse

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

**Immunogen**

A full length human recombinant protein was used as the immunogen for this TL1A antibody.

**Preservatives**

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced) and 0.05% sodium azide

**Storage**

Store the TL1A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

**Note**

For research use only

**Application notes**

Titration of the TL1A antibody may be required for optimal performance.

**Formula**

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide

**Isotype**

Mouse IgG1, kappa

**Clonality**

Monoclonal

**Purity**

Protein G affinity chromatography

**Uniprot ID**

[O95150](#)

**Hazard Information**

This TL1A antibody is available for research use only.

**Entrez**


[9966](#)

**Dilution Range**


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**


12 months from date of receipt.



IHC testing of FFPE human colon carcinom...



IHC testing of FFPE human parathyroid ma...



SDS-PAGE Analysis of Purified, BSA-Free ...