




Product Datasheet

MERTK antibody (orb38645)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	MERTK
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	110249
Uniprot ID	Q12866
NCBI	NP_006334.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:50-10, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of MERTK Antibody ...



Formalin-fixed and paraffin-embedded hum...