

## Product Datasheet

### ZC3H12B antibody (orb386532)

<b>Description</b>	Rabbit polyclonal antibody to ZC3H12B
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	C-terminal region of human ZC3H12B
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Affinity purification
<b>MW</b>	93
<b>Uniprot ID</b>	<a href="#">Q5HYM0</a>
<b>Expiration Date</b>	12 months from date of receipt.