


Product Datasheet

ZNRF3 antibody pair (orb386705)

Description	Rabbit polyclonal antibody to ZNRF3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	N-terminal [ZNRF3 (N) (R2407-1)] or the C-terminal [ZNRF3 (C) (R2407-2)]
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	101
Uniprot ID	Q9ULT6
Dilution Range	WB: 1:1000-:3000, IP: 1:100-200
Expiration Date	12 months from date of receipt.



(A) Western blot analysis of total cell e...