

## Product Datasheet

### His-Tag antibody (orb386759)

<b>Description</b>	Mouse monoclonal antibody to His-Tag
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	His-tag
<b>Preservatives</b>	0.02% Thimerosal
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purification
<b>Expiration Date</b>	12 months from date of receipt.