

Product Datasheet

ZC3H12C antibody (orb386872)

Description	Rabbit polyclonal antibody to ZC3H12C
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal region of human ZC3H12C
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Polyclonal
Purity	Affinity purification
MW	99
Uniprot ID	Q9C0D7
Expiration Date	12 months from date of receipt.