

Product Datasheet

Rabbit IgG Light-Chain antibody (HRP) (orb387000)

Description	Mouse monoclonal antibody against anti-Rabbit IgG Light-Chain conjugated to HRP.
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	rabbit IgG light-chain
Preservatives	0.01% NaN ₃
Storage	-20°C
Note	For research use only
Clonality	Monoclonal
Purity	Affinity purification
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)