

## Product Datasheet

# Rabbit IgG Light-Chain antibody (HRP) (orb387000)

<b>Description</b>	Mouse monoclonal antibody against anti-Rabbit IgG Light-Chain conjugated to HRP.
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC, WB
<b>Immunogen</b>	rabbit IgG light-chain
<b>Preservatives</b>	0.01% NaN <sub>3</sub>
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Affinity purification
<b>Expiration Date</b>	12 months from date of receipt.