

---

## Product Datasheet

### H1 antibody (orb388204)

**Description**

Mouse monoclonal antibody to H1

**Species/Host**

Mouse

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

FC, IF, IHC-P, WB

**Applications**
**Immunogen**

Recombinant full-length human H1F0 Histone H1.0 protein

**Target**

H1

**Preservatives**

200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 1mM PBS with 0.05% rAlbumin &amp; 0.05% azide. Also available WITHOUT rAlbumin &amp; azide at 1.0mg/ml.

**Concentration**

Purified Ab with BSA and Azide at 200ug/ml

**Storage**

Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

**Note**

For research use only

**Application notes**

Histone H1 (Nuclear Marker)

**Isotype**

Mouse / IgG2a, kappa

**Clonality**

Monoclonal

**MW**

~30kDa

**Uniprot ID**


Multiple

**Entrez**
**3005**
**Dilution Range**

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

**Expiration Date**


12 months from date of receipt.



Immunohistochemical staining of rat Panc...



Immunohistochemical staining of human Ov...



Immunohistochemical staining of human To...