
Product Datasheet

HSPD1 antibody (orb388245)

Description

Mouse monoclonal antibody to HSPD1

Species/Host	Mouse
Reactivity	Animal, Bacteria, E. coli, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Plant, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant human HSP60 protein
Target	HSPD1
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody is stable for 24 months. ,Non-hazardous. No MSDS required.
Note	For research use only
Application notes	HSP60 (Heat Shock Protein 60) (Mitochondrial Marker)
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	60kDa
Uniprot ID	P10809
Entrez	3329
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, WB: 0.25-0.5 µg/ml
Expiration Date	12 months from date of receipt.