

Product Datasheet

IGHG antibody (orb388271)

Description

Mouse monoclonal antibody to IGHG

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P

Applications
Immunogen

Purified human IgG

Target

IGHG

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

IgG (Immunoglobulin Gamma Heavy Chain) (B-Cell Marker)

Isotype

Mouse / IgG1, kappa + Mouse / IgG2a, kappa

Clonality

Monoclonal

MW


~75kDa

Uniprot ID
[P01860](#), [P01857](#), [P01859](#), [P01861](#)
Entrez
[3500](#), [3501](#), [3502](#), [3503](#)
Dilution Range

FACS: 1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml

Expiration Date

12 months from date of receipt.


 Immunohistochemical staining of human
To...