


Product Datasheet

IGHG antibody (orb388274)

Description	Rabbit monoclonal antibody to IGHG
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Recombinant full-length human IGHG protein
Target	IGHG
Preservatives	200ug/ml of recombinant MAb purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Form/Appearance	200ug/ml of recombinant MAb purified by Protein A. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	IgG (Immunoglobulin Gamma Heavy Chain) (B-Cell Marker)
Isotype	Rabbit / IgG
Clonality	Recombinant
MW	75kDa
Uniprot ID	P01859, P01861, P01857, P01860
Entrez	3501, 3503, 4, 2, 3500, 3502, 1, 3
Dilution Range	FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human To...