



Product Datasheet

MUC1 antibody (orb388454)

Description	Mouse monoclonal antibody to MUC1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Human milk-fat globule membranes (HMFGM)
Target	MUC1
Preservatives	200ug/ml of Ab purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8 °C. Antibody without azide - store at -20 to -80 °C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	MUC1 / EMA / CD227 (Epithelial Marker)
Isotype	Mouse / IgG2a, lambda
Clonality	Monoclonal
MW	265-400kDa
Uniprot ID	P15941
Entrez	4582
Dilution Range	IHC-P: 0.5-1.0ug/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Ov...



Immunohistochemical staining of human Co...



Immunohistochemical staining of human En...