
Product Datasheet

ST6GAL1 antibody (orb388647)

Description

Mouse monoclonal antibody to ST6GAL1


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P

Immunogen

Nuclei from pokeweed mitogen-stimulated peripheral blood lymphocytes.


Target

ST6GAL1

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by PEG precipitation. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Immunohistochemical staining of Normal h...



Flow cytometric analysis of human Lympho...

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

CDw75 (B-Cell Marker)

Isotype

Mouse / IgM, kappa

Clonality

Monoclonal

MW

46kDa

Uniprot ID
[P15907](#)
Entrez
[6480](#)
Dilution Range

FACS: 0.5-1 µg/million cells, IF/ICC: 0.5-1 µg/ml, IHC-P: 0.5-1.0 µg/ml

Expiration Date

12 months from date of receipt.