


Product Datasheet

HPV-16 antibody (orb389005)

Description	Mouse monoclonal antibody to HPV-16
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Recombinant major capsid protein L1 from human papilloma virus type 16
Target	Not Applicable
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Application notes	HPV-16 (Human Papilloma Virus 16)
Isotype	Mouse / IgG2a, kappa
Clonality	Monoclonal
MW	57kDa
Uniprot ID	Not Applicable
Entrez	Not Applicable
Dilution Range	IF/ICC: 0.5-1 µg/ml, IHC-P (0.5-1.0 µg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Ce...