
Product Datasheet

Hepatocyte Specific Antigen antibody (orb389012)

Description

Mouse monoclonal antibody to Hepatocyte


Species/Host

Mouse

Reactivity

Canine, Human

Conjugation

Unconjugated

Tested Applications

IF, IHC-P



Immunohistochemical staining of human He...

Immunogen

Extract of a formalin-fixed, rejected-allograft of a human liver

Preservatives

200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.

Concentration

Purified Ab with BSA and Azide at 200ug/ml

Storage

Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Application notes

Hepatocyte Specific Antigen (Hepatocellular Marker)

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

Not Known

Uniprot ID

Not Known

Entrez

Not Known

Dilution Range

IF/ICC: 0.5-1.0 µg/ml, IHC-P: 1-2 µg/ml

Expiration Date

12 months from date of receipt.