



Product Datasheet

CXCL10 antibody (orb389312)


Description	CXCL10 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CXCL10. Please contact us for the exact immunogen sequence. The peptide is available as orb13441.
Target	CXCL10
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	11 kDa
Uniprot ID	P02778 , P48973
NCBI	001556
Dilution Range	WB: 1:100-500, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



IHC-P image of human ovarian cancer tiss...



Immunohistochemical staining of paraffin...



Western blot analysis of rat ovary (lane...